



Intermediate capsid protein VP6 Antibody, FITC conjugated

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA321618LC01PGK |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04509 |
| Immunogen | Recombinant Rotavirus A Intermediate capsid protein VP6 protein (1-397AA) |
| Raised In | Rabbit |
| Species Reactivity | Rotavirus A |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Intermediate capsid protein VP6 |
| Immunogen Species | Rotavirus A (strain RVA/Cow/France/RF/1975/G6P6[1]) (RV-A) |
| Research Area | Others |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |