



MMP14 Antibody, Biotin conjugated

Product Code	CSB-PA12389D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50281
Immunogen	Recombinant Human Matrix metalloproteinase-14 protein (112-582AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix metalloproteinase-14 (MMP-14) (EC 3.4.24.80) (MMP-X1) (Membrane-type matrix metalloproteinase 1) (MT-MMP 1) (MTMMP1) (Membrane-type-1 matrix metalloproteinase) (MT1-MMP) (MT1MMP), MMP14
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.