



Toxin zeta Antibody, HRP conjugated

Product Code	CSB-PA363578LB01SMF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A4M2
Immunogen	Recombinant Streptococcus agalactiae Toxin zeta protein (1-287AA)
Raised In	Rabbit
Species Reactivity	Streptococcus agalactiae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin zeta (UDP-N-acetylglucosamine kinase) (UNAG kinase) (EC 2.7.1.176)
Immunogen Species	Streptococcus agalactiae
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.