



RNF7 Antibody, Biotin conjugated

Product Code	CSB-PA019897DD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBF6
Immunogen	Recombinant Human RING-box protein 2 protein (2-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RING-box protein 2 (Rbx2) (CKII beta-binding protein 1) (CKBBP1) (RING finger protein 7) (Regulator of cullins 2) (Sensitive to apoptosis gene protein), RNF7, RBX2 ROC2 SAG
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF7