



Lilrb3 Antibody, Biotin conjugated

Product Code	CSB-PA012941YD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97484
Immunogen	Recombinant Mouse Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (664-841AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3) (Leukocyte immunoglobulin-like receptor 3) (Cell-surface glycoprotein p91) (Paired immunoglobulin-like receptor B) (PIR-B), Lilrb3, Pirb
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Lilrb3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.