



P4hb Antibody, HRP conjugated

Product Code	CSB-PA017342YB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09103
Immunogen	Recombinant Mouse Protein disulfide-isomerase protein (20-509AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase (PDI) (EC 5.3.4.1) (Cellular thyroid hormone-binding protein) (Endoplasmic reticulum resident protein 59) (ER protein 59) (ERp59) (Prolyl 4-hydroxylase subunit beta) (p55), P4hb, Pdia1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	P4hb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.