





NR1I3 Antibody

Product Code	CSB-PA13019A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14994
Immunogen	Recombinant Human Nuclear receptor subfamily 1 group I member 3 protein (2-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 1 group I member 3 (Constitutive activator of retinoid response) (Constitutive active response) (Constitutive androstane receptor) (CAR) (Orphan nuclear receptor MB67), NR1I3, CAR
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NR1I3


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA13019A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA13019A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA13019A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.