



Has2 Antibody

Product Code	CSB-PA010140LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P70312
Immunogen	Recombinant Mouse Hyaluronan synthase 2 protein (67-374AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronan synthase 2 (EC 2.4.1.212) (Hyaluronate synthase 2) (Hyaluronic acid synthase 2) (HA synthase 2), Has2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Has2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.