



RNASEL Antibody, HRP conjugated

Product Code	CSB-PA313223LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05823
Immunogen	Recombinant Human 2-5A-dependent ribonuclease protein (330-741AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-5A-dependent ribonuclease (2-5A-dependent RNase) (EC 3.1.26.-) (Ribonuclease 4) (Ribonuclease L) (RNase L), RNASEL, RNS4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RNASEL