



RNASEL Antibody, HRP conjugated

Product Code	CSB-PA313223LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05823
Immunogen	Recombinant Human 2-5A-dependent ribonuclease protein (330-741AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-5A-dependent ribonuclease (2-5A-dependent RNase) (EC 3.1.26.-) (Ribonuclease 4) (Ribonuclease L) (RNase L), RNASEL, RNS4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RNASEL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.