



Mcl1 Antibody, Biotin conjugated

Product Code	CSB-PA013588HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97287
Immunogen	Recombinant Mouse Induced myeloid leukemia cell differentiation protein Mcl-1 homolog protein (1-308AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Induced myeloid leukemia cell differentiation protein Mcl-1 homolog (Bcl-2-related protein EAT/mcl1), Mcl1
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Mcl1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.