



# cry1Ab Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA363358YC01BDC
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0A370
<b>Immunogen</b>	Recombinant Bacillus thuringiensis subsp. kurstaki Pesticidal crystal protein Cry1Ab protein (1022-1155AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bacillus thuringiensis subsp. kurstaki
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pesticidal crystal protein Cry1Ab (130 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIA(b)), cry1Ab, bt2 cry-1-2 cry1A(b) cryIA(b) cryIC1
<b>Immunogen Species</b>	Bacillus thuringiensis subsp. kurstaki
<b>Research Area</b>	Others
<b>Target Names</b>	cry1Ab
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.