



# TRAC Antibody

Product Code	CSB-PA024144LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01848
Immunogen	Recombinant Human T cell receptor alpha constant protein (1-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T cell receptor alpha constant, TRAC, TCRA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TRAC