




Gnao1 Antibody

Product Code	CSB-PA009593LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18872
Immunogen	Recombinant Mouse Guanine nucleotide-binding protein G(o) subunit alpha protein (2-354AA)
Raised In	Rabbit
Species Reactivity	Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(o) subunit alpha, Gnao1, Gna0 Gnao
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Gnao1

Image



Western Blot

Positive WB detected in: Rat brain tissue, Mouse brain tissue

All lanes: Gnao1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41 kDa

Observed band size: 41 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.