



S100A4 Antibody, Biotin conjugated

Product Code	CSB-PA020632HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26447
Immunogen	Recombinant Human Protein S100-A4 protein (2-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A4 (Calvasculin) (Metastasin) (Placental calcium-binding protein) (Protein Mts1) (S100 calcium-binding protein A4), S100A4, CAPL MTS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	S100A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.