



Serpina1 Antibody, Biotin conjugated

Product Code	CSB-PA021053HD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17475
Immunogen	Recombinant Rat Alpha-1-antiproteinase protein (25-411AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1-antiproteinase (Alpha-1-antitrypsin) (Alpha-1-proteinase inhibitor) (Serpina1), Serpina1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Serpina1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.