



TFF3 Antibody, Biotin conjugated

Product Code	CSB-PA023433YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07654
Immunogen	Recombinant Human Trefoil factor 3 protein (22-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trefoil factor 3 (Intestinal trefoil factor) (hITF) (Polypeptide P1.B) (hP1.B), TFF3, ITF TFI
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TFF3