



IL1RL1 Antibody

Product Code	CSB-PA07209A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01638
Immunogen	Recombinant Human Interleukin-1 receptor-like 1 protein (19-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 receptor-like 1 (Protein ST2), IL1RL1, DER4 ST2 T1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IL1RL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.