



Parvalbumin beta Antibody, Biotin conjugated

Product Code	CSB-PA365766HD01GAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02622
Immunogen	Recombinant Gadus morhua subsp Parvalbumin beta protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Gadus morhua subsp
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Parvalbumin beta (Allergen Gad c I) (Allergen M) (allergen Gad c 1)
Immunogen Species	Gadus morhua subsp. callarias (Baltic cod) (Gadus callarias)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.