



COL12A1 Antibody

Product Code	CSB-PA005719LA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13944
Immunogen	Recombinant Chicken Collagen alpha-1(XII) chain protein (139-311AA)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-1(XII) chain (Fibrochimerin), COL12A1
Immunogen Species	Chicken
Research Area	Others
Target Names	COL12A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.