



Ifne Antibody, HRP conjugated

Product Code	CSB-PA800374LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q80ZF2
Immunogen	Recombinant Mouse Interferon epsilon protein (22-192AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon epsilon (IFN-epsilon) (Interferon epsilon-1) (Interferon tau-1), Ifne, Ifne1 Ifnt1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Ifne
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.