



TNFAIP6 Antibody, Biotin conjugated

Product Code	CSB-PA023959LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98066
Immunogen	Recombinant Human Tumor necrosis factor-inducible gene 6 protein (18-277AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor-inducible gene 6 protein (Hyaluronate-binding protein) (TNF-stimulated gene 6 protein) (TSG-6) (Tumor necrosis factor alpha-induced protein 6) (TNF alpha-induced protein 6), TNFAIP6, TSG6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFAIP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.