



# TNFAIP6 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA023959LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P98066
<b>Immunogen</b>	Recombinant Human Tumor necrosis factor-inducible gene 6 protein (18-277AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tumor necrosis factor-inducible gene 6 protein (Hyaluronate-binding protein) (TNF-stimulated gene 6 protein) (TSG-6) (Tumor necrosis factor alpha-induced protein 6) (TNF alpha-induced protein 6), TNFAIP6, TSG6
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Immunology
<b>Target Names</b>	TNFAIP6
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.