



Genome polyprotein Antibody

Product Code	CSB-PA714772LA01HYH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66478
Immunogen	Recombinant Human enterovirus 71 Genome polyprotein protein (566-862AA)
Raised In	Rabbit
Species Reactivity	Human enterovirus 71
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Genome polyprotein [Cleaved into: P3; Protein 3AB; P2; P1; Capsid protein VP0 (VP4-VP2); Capsid protein VP4 (P1A) (Virion protein 4); Capsid protein VP2 (P1B) (Virion protein 2); Capsid protein VP3 (P1C) (Virion protein 3); Capsid protein VP1 (P1D) (Virion protein 1); Protease 2A (P2A) (EC 3.4.22.29) (Picornain 2A) (Protein 2A); Protein 2B (P2B); Protein 2C (P2C) (EC 3.6.1.15); Protein 3A (P3A); Viral protein genome-linked (VPg) (Protein 3B) (P3B); Protein 3CD (EC 3.4.22.28); Protease 3C (P3C) (EC 3.4.22.28); RNA-directed RNA polymerase (RdRp) (EC 2.7.7.48) (3D polymerase) (3Dpol) (Protein 3D) (3D)]
Immunogen Species	Human enterovirus 71
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.