



ACHE Antibody

| | |
|---------------------|--|
| Product Code | CSB-PA001154EA01BO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23795 |
| Immunogen | Recombinant Bovine Acetylcholinesterase protein (31-613AA) |
| Raised In | Rabbit |
| Species Reactivity | Bovine |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Acetylcholinesterase (AChE) (EC 3.1.1.7), ACHE |
| Immunogen Species | Bos taurus (Bovine) |
| Research Area | Neuroscience |
| Target Names | ACHE |