



ACHE Antibody

Product Code	CSB-PA001154EA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23795
Immunogen	Recombinant Bovine Acetylcholinesterase protein (31-613AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetylcholinesterase (AChE) (EC 3.1.1.7), ACHE
Immunogen Species	Bos taurus (Bovine)
Research Area	Neuroscience
Target Names	ACHE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.