




S100A6 Antibody

Product Code	CSB-PA020634HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06703
Immunogen	Recombinant Human Protein S100-A6 protein (1-90AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A6 (Calcyclin) (Growth factor-inducible protein 2A9) (MLN 4) (Prolactin receptor-associated protein) (PRA) (S100 calcium-binding protein A6), S100A6, CACY
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	S100A6

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA020634HA01HU at dilution of 1:100