



# Minor capsid protein VP2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA355948LC01JAK
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P03095
<b>Immunogen</b>	Recombinant JC polyomavirus Minor capsid protein VP2 protein (2-344AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	JC polyomavirus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Minor capsid protein VP2 (Minor structural protein VP2)
<b>Immunogen Species</b>	JC polyomavirus (JCPyV) (JCV)
<b>Research Area</b>	Others