



Minor capsid protein VP2 Antibody, HRP conjugated

Product Code	CSB-PA355948LB01JAK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03095
Immunogen	Recombinant JC polyomavirus Minor capsid protein VP2 protein (2-344AA)
Raised In	Rabbit
Species Reactivity	JC polyomavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Minor capsid protein VP2 (Minor structural protein VP2)
Immunogen Species	JC polyomavirus (JCPyV) (JCV)
Research Area	Others