



Minor capsid protein VP2 Antibody, HRP conjugated

Product Code	CSB-PA355948LB01JAK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03095
Immunogen	Recombinant JC polyomavirus Minor capsid protein VP2 protein (2-344AA)
Raised In	Rabbit
Species Reactivity	JC polyomavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Minor capsid protein VP2 (Minor structural protein VP2)
Immunogen Species	JC polyomavirus (JCPyV) (JCV)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.