




OGDH Antibody

Product Code	CSB-PA016307LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02218
Immunogen	Recombinant Human 2-oxoglutarate dehydrogenase, mitochondrial protein (41-423AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-oxoglutarate dehydrogenase, mitochondrial (EC 1.2.4.2) (2-oxoglutarate dehydrogenase complex component E1) (OGDC-E1) (Alpha-ketoglutarate dehydrogenase), OGDH
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OGDH

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA016307LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.