



LGALS9 Antibody, HRP conjugated

Product Code	CSB-PA02649B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00182
Immunogen	Recombinant Human Galectin-9 protein (1-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-9 (Gal-9) (Ecalectin) (Tumor antigen HOM-HD-21), LGALS9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LGALS9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.