



E2 Antibody, FITC conjugated

Product Code	CSB-PA361949LC01HMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06790
Immunogen	Recombinant Human papillomavirus type 18 Regulatory protein E2 protein (1-365AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 18
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulatory protein E2, E2
Immunogen Species	Human papillomavirus type 18
Research Area	Others
Target Names	E2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.