



## E2 Antibody

Product Code	CSB-PA361949LA01HMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06790
Immunogen	Recombinant Human papillomavirus type 18 Regulatory protein E2 protein (1-365AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 18
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulatory protein E2, E2
Immunogen Species	Human papillomavirus type 18
Research Area	Others
Target Names	E2