



# S100g Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA020644LB01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P02634
<b>Immunogen</b>	Recombinant Rat Protein S100-G protein (2-79AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein S100-G (9 kDa CaBP) (Calbindin-D9k) (Cholecalciferol) (S100 calcium-binding protein G) (Vitamin D-dependent calcium-binding protein, intestinal) (CABP), S100g, Calb3 S100d
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Others
<b>Target Names</b>	S100g
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.