



S100g Antibody, HRP conjugated

Product Code	CSB-PA020644LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02634
Immunogen	Recombinant Rat Protein S100-G protein (2-79AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-G (9 kDa CaBP) (Calbindin-D9k) (Cholecalciferol-binding protein G) (Vitamin D-dependent calcium-binding protein, intestinal) (CABP), S100g, Calb3 S100d
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	S100g