



MMP3 Antibody, Biotin conjugated

Product Code	CSB-PA07449D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08254
Immunogen	Recombinant Human Stromelysin-1 protein (100-477AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stromelysin-1 (SL-1) (EC 3.4.24.17) (Matrix metalloproteinase-3) (MMP-3) (Transin-1), MMP3, STMY1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.