




AAO Antibody

Product Code	CSB-PA333022LA01DZK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24792
Immunogen	Recombinant Cucurbita maxima L-ascorbate oxidase protein (31-579AA)
Raised In	Rabbit
Species Reactivity	Cucurbita maxima
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-ascorbate oxidase (ASO) (Ascorbase) (EC 1.10.3.3), AAO
Immunogen Species	Cucurbita maxima
Research Area	Others
Target Names	AAO

Image



Western Blot

Positive WB detected in: Cucurbita maxima (40µg, 20µg, 10µg, 5µg)

All lanes: AAO antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 65 kDa

Observed band size: 65 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.