



Gpc4 Antibody

Product Code	CSB-PA009706EA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51655
Immunogen	Recombinant Mouse Glypican-4 protein (530-557AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glypican-4 (K-glypican) [Cleaved into: Secreted glypican-4], Gpc4
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Gpc4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.