



# KLK1 Antibody, HRP conjugated

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA012446LB01HU                                                                                |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| Uniprot No.         | P06870                                                                                            |
| Immunogen           | Recombinant Human Kallikrein-1 protein (25-262AA)                                                 |
| Raised In           | Rabbit                                                                                            |
| Species Reactivity  | Human                                                                                             |
| Tested Applications | ELISA                                                                                             |
| Form                | Liquid                                                                                            |
| Conjugate           | HRP                                                                                               |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                  |
| Purification Method | >95%, Protein G purified                                                                          |
| Isotype             | IgG                                                                                               |
| Clonality           | Polyclonal                                                                                        |
| Alias               | Kallikrein-1 (EC 3.4.21.35) (Kidney/pancreas/salivary gland kallikrein) (Tissue kallikrein), KLK1 |
| Immunogen Species   | Homo sapiens (Human)                                                                              |
| Research Area       | Cardiovascular                                                                                    |
| Target Names        | KLK1                                                                                              |