



KLK1 Antibody, Biotin conjugated

Product Code	CSB-PA012446LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06870
Immunogen	Recombinant Human Kallikrein-1 protein (25-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-1 (EC 3.4.21.35) (Kidney/pancreas/salivary gland kallikrein) (Tissue kallikrein), KLK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	KLK1