



APOA1 Antibody, HRP conjugated

Product Code	CSB-PA001913LB01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15497
Immunogen	Recombinant Bovine Apolipoprotein A-I protein (25-265AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein A-I (Apo-AI) (ApoA-I) (Apolipoprotein A1) [Cleaved into: Proapolipoprotein A-I (ProapoA-I); Truncated apolipoprotein A-I], APOA1
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	APOA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.