



E7 Antibody

Product Code	CSB-PA361948LA01HMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06788
Immunogen	Recombinant Human papillomavirus type 18 Protein E7 protein (1-105AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 18
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E7, E7
Immunogen Species	Human papillomavirus type 18
Research Area	Others
Target Names	E7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.