



Angptl8 Antibody

Product Code	CSB-PA844436LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R1L8
Immunogen	Recombinant Mouse Betatrophin protein (16-198AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-like protein 8 (Betatrophin) (Lipasin) (Refeeding-induced fat and liver protein), Angptl8, Gm6484 R1f1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Angptl8