



Angptl8 Antibody, Biotin conjugated

Product Code	CSB-PA844436LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R1L8
Immunogen	Recombinant Mouse Betatrophin protein (16-198AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-like protein 8 (Betatrophin) (Lipasin) (Refeeding-induced fat and liver protein), Angptl8, Gm6484 R1f1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Angptl8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.