




Afp Antibody

Product Code	CSB-PA001421LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02772
Immunogen	Recombinant Mouse Alpha-fetoprotein protein (18-605AA)
Raised In	Rabbit
Species Reactivity	Mouse, Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin), Afp
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Afp

Image



Western blot
 All lanes: Afp antibody at 10µg/ml
 Lane 1: HepG2 whole cell lysate
 Lane 2: Rat liver tissue
 Lane 3: 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 68 kDa
 Observed band size: 68 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.