



CXCL8 Antibody, HRP conjugated

Product Code	CSB-PA15789B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08317
Immunogen	Recombinant Gallus gallus Interleukin-8 protein (17-102AA)
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-8 (IL-8) (9E3) (C-X-C motif chemokine 8) (CEF-4) (Chemokine (C-X-C motif) ligand 8) (Embryo fibroblast protein 1) (EMF-1), CXCL8, CEF4 EMF1 IL8
Immunogen Species	Gallus gallus (Chicken)
Research Area	Others
Target Names	CXCL8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.