



AQP4 Antibody

Product Code	CSB-PA001964GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	For research use only. Not for human, diagnostic or therapeutic use.
Uniprot No.	P55087
Immunogen	Human AQP4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	HMIWC2; MGC22454; MIWC; aquaporin 4
Immunogen Species	Homo sapiens (Human)
Target Names	AQP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.