





PNMA5 Antibody

Product Code	CSB-PA853483HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PV4
Immunogen	Recombinant Human Paraneoplastic antigen-like protein 5 protein (4-442AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paraneoplastic antigen-like protein 5 (Tumor antigen BJ-HCC-25), PNMA5, KIAA1934
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PNMA5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA853483HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA853483HA01HU at dilution of 1:100