



DUSP13 Antibody, FITC conjugated

Product Code	CSB-PA734923HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6B8I1
Immunogen	Recombinant Human Dual specificity protein phosphatase 13 isoform A protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 13 isoform A (DUSP13A) (EC 3.1.3.16) (EC 3.1.3.48) (Branching-enzyme interacting DSP) (Muscle-restricted DSP) (MDSP), DUSP13, BEDP DUSP13A MDSP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DUSP13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.